



# DATA SHEET

ADJUSTABLE CIRCULAR HAT. REF: LQ1016

**DESCRIPTION:**

Cap made of non-woven, non-sterile polypropylene, which acts as a protective barrier to prevent microorganisms that are shed from the hair. can cause contamination problems.

**FEATURES:**

17 gram SSnon woven fabric with 0.5 cm elastic

**TECHNICAL SPECIFICATIONS:**

28 cm diameter with 0.5 cm elastic

**HOW TO USE:**

Place on the head making sure to completely cover the hair (including the equillo). This hat should be put on before dressing, to protect the clothing from the contamination by hair.

**COMMERCIAL PRESENTATION:**

12 units per package

**STORAGE:**

This product should be stored in a clean, cool and contamination-free place. Temperature 5 to 30 °C

**WARRANTY:**

15 días por defectos de fabricación

**FINAL DISPOSITION:**

Disposable product for single use. Do not reuse. It must be disposed of in containers where medical devices are stored.

**APPLICATIONS:**

- Reduces contamination of the surgical field by microbes shed from the hair and scalp.
- Prevents the loss of desquamated cells from the scalp that may contain bacteria, and be a source of contamination for the patient or the surfaces and instruments in the operating room.
- They are used mainly in operating rooms, in units with critically ill patients and in strict isolation.

**CARE AND CLEANING:**

- Wash your hands prior to contact with the garment.
- Put on the cap covering all the hair, in case of having long hair, it must be collected previously for better support.
- Discard after use.
- The materials used do not cause any allergic reaction, are not toxic or cytotoxic and do not lead to sensitization or irritation.